

IPC-09 DEAD WEIGHT VACUUM GAUGE TESTER



RANGE OF EXPERIMENTS TO BE CARRIED OUT:

1. To calibrate the given Vacuum gauge

TECHNICAL DESCRIPTION :

The setup consist of screw pump to generate the vacuum. The free piston assembly is connected to the screw pump. The set of dead weights is attached to the piston assembly. The whole setup is mounted on sturdy MS plate with good color combination.

DIMENSIONS AND WEIGHT :

Size :0.5 m.(L)x 0.5m(W) X 0.4m (H)

Weight :Approx. 3 Kg

SERVICE REQUIRED :

SCOPE OF DELIVERY:

1. Experimental Setup
2. Instructional Manual

1, Camp's Corner, Nr. Narhari Hospital, Fatehgunj, Vadodara - 390 002, Gujarat, India.

Tell No. +91 265 750186, Cell No. +91 9727759429

info@fadaklabequipments.com