

SHM-04 BATCH CRYSTALIZER



RANGE OF EXPERIMENTS TO BE CARRIED OUT:

- 1. To study the crystallization Process
- 2. To calculate the percentage yield of crystallization

TECHNICAL DESCRIPTION :

- i. It consists of vessel SS-316 of 1 Ltr capacity provided with bottom discharge valve.
- ii. Water Bath of 5 ltr capacity
- iii. Agitator: SS-316 rotated by fraction horse power, single phase motor.
- iv. Cooling water arrangement : BY REFRIGERATOIN SYSTEM

DIMENSIONS AND WEIGHT :

Size :1.5 m.(L)x 1m(W) X 1.5m (H) Weight :Approx. 50 Kg

SERVICE REQUIRED : 230 v Ac Supply 50 Hz

SCOPE OF DELIVERY:

- 1. Experimental Setup
- 2. Instructional Manual

1, Camp's Corner, Nr. Narhari Hospital, Fatehgunj, Vadodara - 390 002, Gujarat, India. Tell No. +91 265 750186, Cell No. +91 9727759429 info@fadaklabequipments.com