

IPC-004 TIME CONSTANT OF MANOMETER



RANGE OF EXPERIMENTS TO BE CARRIED OUT:

1. To Determine the time constant of given Manometer

TECHNICAL DESCRIPTION :

Here Three Different lengths of manometer are selected & used to measure the manometer difference. An air filter regulator is used to provide air at constant pressure.

DIMENSIONS AND WEIGHT :

Size :0.7 m.(L)x 0.5m(W) X 1.2m (H)

Weight :Approx. 20 Kg

SERVICE REQUIRED :

Air Supply : @2 kg/cm²