

FAL-035 STOKES LAW APPARATUS



RANGE OF EXPERIMENTS TO BE CARRIED OUT:

1. To Verify the stokes law

TECHNICAL DESCRIPTION :

The setup consist of Consist of glass tubes of suitable length and diameter. This glass tubes are filled with oil. Small balls of different diameters are provided.The whole setup mounted on suitable stand .

DIMENSIONS AND WEIGHT :

Size :0.6 m.(L)x 0.5 m(W) X 1.0m (H)

Weight :Approx. 20Kg

SERVICE REQUIRED :

SCOPE OF DELIVERY:

1. Experimental Setup
2. Instructional Manual

1, Camp's Corner, Nr. Narhari Hospital, Fatehgunj, Vadodara - 390 002, Gujarat, India.

Tell No. +91 265 750186, Cell No. +91 9727759429

info@fadaklabequipments.com